

For European Union Countries

Type \ Specs.	Model	Photo	Frequency (MHz)	Requied Amplifier	Sterilization method	Connecting Dopplers				
						ES-100VX	ES-100V3	Bidop 7	Smartdop 45	Smartdop XT
Vascular probe: Standard	BT*M05S8C* [†]		4/ 5/ 8/ 10	-	-	✓	✓	-	✓	4/ 5/ 8
	BT*M05S8A* [†]									
	ST*M05S8C* [†]		4/ 5/ 8/ 10	-	-	-	-	✓	-	-
	ST*M05S8A* [†]									
2MHz probe: Standard	BT2M20S8C* [†]		2	-	-	✓	✓	-	✓	-
	BT2M20S8A* [†]									
PPG probe	PG-21		-	-	-	-	✓	-	✓	-
	PG-30		-	-	-	-	-	✓	-	-

*[†]: * The last word “C” represents Curl type cable and “A” represents Straight type one.

*Frequency